


Face Shield RFE078-1159

Made in the UK

Having dedicated our efforts to helping the NHS protect its front-line staff, we are now able to provide PPE to organisations looking to ensure the safety of customer facing workers



From a leading provider of PPE to the NHS throughout the Covid-19 Pandemic, this Face Shield is suitable for medical and non-medical environments

- Tested to the COVID-19 requirements of BS EN166:2001
- The frameless screen, manufactured from lightweight, optically clear PET, provides a distortion free, unobstructed field of vision
- Screen size: 320 x 230mm
- Anti-fog and anti-static treatments ensure vision remains clear
- The foam forehead rest ensures comfort even during extended wearing
- The size adjustable headband provides a comfortable fit for all head sizes.
- There is room behind the visor for the wearing of a face mask and spectacles
- Branded versions of the shield are available

Minimum Order Quantity: 5000

Price: £1.30

For more information:

Adam Lyon

Tel: 07895 200843

e-mail: adam.lyon@rociallehealthcare.com

Sarah Staufenbiel

Tel: 07458 086698

e-mail: sarah.staufenbiel@rociallehealthcare.com